

HPLC Vacuum Degasser Module Contamination (sent to us for testing & evaluation):



*Many of these parts were purchased or removed from degassers sold by sellers on the Internet, described as refurbished or “like new”. To find out if your degasser is contaminated or to prevent contamination of the HPLC flow path, please have your HPLC degasser professionally serviced by Chiralizer Services, LLC at least every 5 years [More info: www.chiralizer.com ; email sales@chiralizer.com or sales@hplctools.com].

Copyright 2020, Chiralizer Services, LLC (USA). All Rights Reserved.